







## THE I.S.S. RANGE OF FILTER PIPETTE TIPS

Rack Packed Filter Tips				
Code	Volume $\mu$ l	Colour	Case Qty	Rack Qty
FT0015	0.5-8	White	960	96
FT0020	1-20	White	960	96
FT0050	1-100	White	960	96
FT0150	1-140	White	960	96
FT0250	1-230	White	960	96
FT1000	200-1000	White	960	96

FT0015	FT0020	FT0050	FT00150	FT0250	FT1000
					
0.5-10 microlitres	1-30 microlitres	1-100 microlitres	1-180 microlitres	1-230 microlitres	200-1000 microlitres

## COMPATIBILITY CHART

Eppendorf * Research		✓	✓		✓	✓
Eppendorf * Reference		✓	✓		✓	✓
Gilson * Pipetman P2	✓					
Gilson * Pipetman P10	✓					
Gilson * Pipetman P20		✓				
Gilson * Pipetman P100			✓			
Gilson * Pipetman P200				✓	✓	
Gilson * Pipetman P1000						✓
Biohit *		✓	✓		✓	✓
Finnpipettes * Digital		✓	✓		✓	✓

\* ®Registered Trademarks