

Orange Screw Cap Paediatric Blood Collection Tube (0.5ml Lithium Heparin)(Pack 1000)

Short Description: Paediatric blood collection tube orange screw cap with lithium heparin 0.5ml filled.

SKU: PAED95

Case Quantity: 1000

Description: Blood collection tube filled with 0.5ml lithium heparin, intended for blood samples. Lithium heparin anticoagulant are commonly used for chemical measurements and clinical biochemistry tests. Due to its low cation concentration, minimal chelating properties and minimal interference with water, heparin is an ideal anticoagulant for chemical analysis in plasma and blood tests.

The orange screw cap tube will allow 0.5ml volume of blood to be collected. This long narrow tube is ideal for sample collection and processing. The tube is CE marked in accordance with EU directive 98/79/EC and is tested to ISO 6710:2017, it is a single-use containers for human venous blood specimen collection. It is specifically designed for paediatric use and colour coded for standard compliance and ease of recognition.

Product Features:

- Tube Type: Paediatric Blood Collection Tube
- Anticoagulant: Lithium Heparin (ISO 6710:2017)
- Sample Volume: 0.5ml
- Cap Type: Screw Cap
- Cap Colour: Orange
- Pierceable: No
- Base: Skirted Conical Base
- Product Material: Polypropylene
- Product Colour: Transparent
- Label Type: Paper

Key Features:

TUBE TYPE: Paediatric Blood Collection Tube

ANTICOAGULANT: 0.5ml Lithium Heparin

CAP COLOUR: Orange

CAP TYPE: Screw Cap

SKU: PAED95

UNIT SIZE: 1000

Website URL Link

<https://www.intscientific.com/product/orange-screw-cap-paediatric-blood-collection-tube-0-5ml-lithium-heparin/>

Email: sales@intscientific.com

W: www.intscientific.com

T: +44(0)1274 720070