


Bags and Bins



Bags / Self Seal

Minigrip® Bags, Plain & Write On

- High quality resealable bags with or without write on surface.



Plain Minigrip® Bag	Write on Minigrip® Bag	Dimensions mm	Case Quantity
MGB001	MGB014	40 x 65	1000
MGB002	MGB015	60 x 60	1000
MGB003	MGB016	75 x 85	1000
MGB004	MGB017	90 x 115	1000
MGB005	MGB018	115 x 115	1000
MGB006	MGB019	100 x 140	1000
MGB007	MGB020	140 x 140	1000
MGB008	MGB021	75 x 190	1000
MGB009	MGB022	125 x 190	1000
MGB010	MGB023	190 x 190	1000
MGB011	MGB024	150 x 230	1000
MGB012	MGB025	205 x 280	1000
MGB036	MGB040	230 x 325	1000
MGB044	MGB051	255 x 355	1000
MGB045	-	280 x 405	1000

Bags / Self Seal with Document Pouch


- Specimen bags with separate document pouch with or without biohazard marking.



Code	Description	Bag Dimensions mm	Pouch Dimensions mm	Case Quantity
MGB026	Specimen Bag Plain	125 x 140	200	1000
MGB038	Specimen Bag Biohazard	125 x 140	200	1000
MGB039	Specimen Bag Plain	140 x 140	200	1000
MGB027	Specimen Bag Biohazard	140 x 140	200	1000
MGB028	Specimen Bag Plain	150 x 140	200	1000
MGB029	Specimen Bag Biohazard	150 x 140	200	1000
MGB050	Specimen Bag Biohazard	250 x 300	250	1000
MGB052	Specimen Bag Biohazard	150 x 254	242	500

Autoclave Bags


- High temperature resistance.
- Bags carry hazard warnings.
- Manufactured from heavy gauge polypropylene for superior performance.
- Autoclavable to 145 degrees centigrade.



Code	Description	Dimensions mm	Case Quantity
ACB001	Autoclave Bag. Small	310 x 662	200
ACB003	Autoclave Bag. Medium	413 x 596	200
ACB005	Autoclave Bag. Large	608 x 816	200

ACT001	Autoclave Tape. 50 metre roll	19mm x 50m	1
--------	-------------------------------	------------	---

Homogeniser Bags



- Sterile homogenizer bags for use with the models listed below.

Code	Description	Dimensions mm	Case Quantity
HMB005	Plain Bag for 80ml Model. Inner Sleeves of 50	152 x 101	2000
HMB003	Plain Bag for 400ml Model. Inner Sleeves of 50	305 x 184	500
HMB004	Curved Bag for 400ml Model. Inner Sleeves of 50	305 x 184	500
HMB007	Twist Seal Bag for 400ml Model. Inner Sleeves of 10	305 x 184	250
HMB008	Strainer Bag for 400ml Model. Inner Sleeves of 10	305 x 184	200
HMB009	Plain Bag for 3500ml Model. Inner Sleeves of 50	515 x 380	250
SS0001	Sampling Straw	175 x 4.75	2000

Sweetie Jars

- Blow moulded, clear sweetie jars.
- For the containment of soft and non-liquid waste prior to safe disposal.



Code	Description	Capacity	Quantity
SWE003	Small Sweetie Jar	2.5 litres	1
SWE002	Large Sweetie Jar	4.5 litres	1

Sharps Bin

- Please call our office for technical specification and extended range.

Code	Description	Quantity
BIN001	4 litre Sharps Bin	1



Warning Tapes

- Self adhesive tape sold in 66m rolls, pre-printed with either 'Danger Of Infection' or 'Pathological Specimen - Fragile Handle With Care'.



Code	Description	Dimensions	Case Quantity
TAPE01	Pathological Specimen - Fragile Handle With Care Tape	25mm x 66m	1 Roll
TAPE02	Danger of Infection Tape. Black on Yellow plus symbol	25mm x 66m	1 Roll