

Blood Tubes



Plastic Screw Cap Blood Collection Tubes with Anticoagulants (5-10ml)

Pre-Filled With Anticoagulant

- A range of screw capped 5 & 10ml blood collection tubes available with or without labels.
- All tubes are fitted with rigid screw caps for leak free specimens.
- Available in polystyrene or clarified polypropylene.
- CE marked in accordance with EU directive 98/79/EC.
- Tested to ISO 6710:2017.



* Contact us if you do not see the blood tube you require and we will endeavour to make it for you.

Code	Volume ml	Anticoagulant	Colour Code	Material	Case Quantity
PBT064	2.5	Di Potassium EDTA	● Pink	Polystyrene	1000
PBT065	2.5	Di Potassium EDTA	● Pink	Polypropylene	1000
PBT068	4	Di Potassium EDTA	● Pink	Polystyrene	1000
PBT089	4	Di Potassium EDTA	● Pink	Polypropylene	1000
PBT088	10	Di Potassium EDTA	● Pink	Polypropylene	500
PBT108	5	Lithium Heparin	● Orange	Polystyrene	1000
PBT107	5	Lithium Heparin	● Orange	Polypropylene	1000
PBT080	10	Lithium Heparin	● Orange	Polystyrene	500
PBT110	10	Lithium Heparin	● Orange	Polypropylene	500
PBT111	10	Lithium Heparin & Beads	● Orange	Polypropylene	500
PBT042	2.5	Sodium Citrate 3.2% (9NC) (Dry).	● Purple	Polystyrene	1000
PBT097	2.5	Fluoride / Oxalate	● Yellow	Polystyrene	1000
PBT096	2.5	Fluoride / Oxalate	● Yellow	Polypropylene	1000

Plastic Screw Cap Blood Collection Tubes without Anticoagulants (5-10ml)

Empty Tubes

- A range of screw capped 5 & 10ml blood collection tubes available with or without labels.
- All tubes are fitted with rigid screw caps for leak free specimens.
- Gamma irradiated (>25Kgy) or non-sterile
- CE marked in accordance with EU directive 98/79/EC.
- Leak tested to ISO 6710:2017.



Code	Volume ml	Anticoagulant	Colour Code	Label	Material	Gamma Irradiated	Case Qty
PBT081	5	Empty	○ White	✓	Polystyrene		1000
PBT128	5	Empty	○ White	✓	Polypropylene		1000
PBT199	5	Empty	○ White	✓	Polystyrene	✓	1000
PBT200	5	Empty	○ White	✓	Polypropylene	✓	1000
PBT166	5	Empty	○ White		Polystyrene		1000
PBT167	5	Empty	○ White		Polypropylene		1000
PBT084	5	Empty	○ White		Polystyrene	✓	1000
PBT201	5	Empty	○ White		Polypropylene	✓	1000
PBT082	10	Empty	○ White	✓	Polystyrene		500
PBT129	10	Empty	○ White	✓	Polypropylene		500
PBT202	10	Empty	○ White	✓	Polystyrene	✓	500
PBT203	10	Empty	○ White	✓	Polypropylene	✓	500
PBT168	10	Empty	○ White		Polystyrene		500
PBT169	10	Empty	○ White		Polypropylene		500
PBT204	10	Empty	○ White		Polystyrene	✓	500
PBT205	10	Empty	○ White		Polypropylene	✓	500

With Clotting Accelerator

- Tubes as above but treated with clotting accelerator.
- Also available with serum separation gel or beads.
- Leak tested to ISO 6710:2017.



Code	Volume ml	Additive	Colour Code	Label	Material	Case Qty
PBT134	5	Treated for Clotted Blood	○ White	✓	Polypropylene	1000
PBT130	5	Treated with Beads	○ White	✓	Polypropylene	1000
PBT145	5	Treated with Separation Gel	● Green	✓	Polypropylene	300
PBT133	10	Treated for Clotted Blood	○ White	✓	Polypropylene	500
PBT131	10	Treated with Beads	○ White	✓	Polypropylene	500
PBT207	10	Treated with Separation Gel	● Green	✓	Polypropylene	500

Glass Narrow & Wide Neck Tubes (5-10ml)

Without Anticoagulant

- A range of screw capped 5 & 10ml blood collection tubes available with or without labels.
- Gamma irradiated (at >25Kgy) or non sterile.
- Fitted with rigid screw caps for leak free specimen collection.
- Available in narrow or wide neck styles.
- 10ml tubes are manufactured in a special process which prevents clot retraction problems.
- CE marked in accordance with EU directive 98/79/EC.
- Leak tested to ISO 6710:2017.



* Contact us if you do not see the blood tube you require and we will endeavour to make it for you.

Code	Volume ml	Description/Neck Style	Label	Gamma Irradiated	Case Quantity
GBT016	5	Narrow			1000
GBT045	5	Narrow	✓		1000
GBT047	5	Narrow	✓	✓	1000
GBT049	5	Narrow		✓	1000
GBT029	5	Wide			1000
GBT028	5	Wide	✓		1000
GBT030	5	Wide	✓	✓	1000
GBT031	5	Wide		✓	1000
GBT021	10	Narrow			500
GBT046	10	Narrow	✓		500
GBT048	10	Narrow	✓	✓	500
GBT050	10	Narrow		✓	500
GBT036	10	Wide			500
GBT035	10	Wide	✓		500
GBT033	10	Wide	✓	✓	500
GBT032	10	Wide		✓	500

Screw Cap Paediatric Blood Collection Tubes (0.5-1.3ml)

Standard Cap

- Screw cap tubes for 0.5 - 1.3ml of blood.
- These long narrow tubes aid storage and "harvest".
- We are able to provide these tubes with a range of additives or treatments.
- CE marked in accordance with EU directive 98/79/EC.
- Tested to ISO 6710:2017.
- Contact us if you do not see the tube you require and we will endeavour to make it for you.



We can manufacture these tubes in alternative colour codes.

Code	Volume ml	Anticoagulant	Colour Code	Pierceable	Material	Case Qty
PAED95	0.5	Lithium Heparin	● Orange	N	Polypropylene	1000
PAED06	1.3	Lithium Heparin	● Orange	N	Polypropylene	1000
PAED20	1.3	Lithium Heparin with Separation Gel	● Orange	N	Polypropylene	1000
PAED01B	0.5	Di Potassium EDTA	● Pink	N	Polypropylene	1000
PAED77	0.8	Di Potassium EDTA	● Pink	N	Polypropylene	1000
PAED02B	1.3	Di Potassium EDTA	● Pink	N	Polypropylene	1000
PAED09	1.3	Treated for Clotted Blood	○ White	N	Polypropylene	1000
PAED04	1.3	Fluoride/Oxalate	● Yellow	N	Polypropylene	1000
PAED76	1.3	Fluoride/Potassium EDTA	● Grey	N	Polypropylene	1000
PAED08	1.3	Plain 'Z'	○ White	N	Polypropylene	1000
PAED19	1.3	Plain, 'Z', Irradiated	○ White	N	Polypropylene	1000
PAED21	1.3	Separation Gel for Clotted Blood	● Brown/White	N	Polypropylene	1000
PAED30	1.3	Sodium Citrate 3.2% (9NC).	● Mauve	N	Polypropylene	100

Pierceable Blood Collection Tubes (0.5-1.3ml)

- Screw caps feature a pierceable membrane.
- Carrier tubes with special 'click fit' feature convert the product to a standard diameter and height for use with automated sampling equipment.
- CE marked in accordance with EU directive 98/79/EC.
- Tested to ISO 6710:2017.
- Contact us if you do not see the tube you require and we will endeavour to make it for you.



We can manufacture these tubes in alternative colour codes.

Code	Volume ml	Anticoagulant	Colour Code	Pierceable	Material	Case Quantity
PAED44	0.8	Sodium Citrate 3.2% (9NC).	● Mauve	Y	Polypropylene	100
PAED71	1.3	Lithium Heparin	● Orange	Y	Polypropylene	1000
PAED49	0.8	Di Potassium EDTA	● Pink	Y	Polypropylene	1000
PAED43	1.3	Di Potassium EDTA	● Pink	Y	Polypropylene	1000
PTT033	5.5	Carrier Tube			Polypropylene	2000

Push Cap Paediatric Blood Collection Tubes (0.5-2.0ml)

- Push cap tubes for 0.5 - 2.0ml of blood.
- Moulded in crystal clear polystyrene.
- CE marked in accordance with EU directive 98/79/EC.
- Tested to ISO 6710:2017.
- We are able to provide these tubes with a range of additives or treatments.



Code	Volume ml	Anticoagulant	Colour Code	Material	Case Qty
PBT0591	1	Empty	○ White	Polystyrene	1000
PBT059	2	Empty	○ White	Polystyrene	1000
PBT062	0.5	Di Potassium EDTA	● Pink	Polystyrene	1000
PBT054	1	Di Potassium EDTA	● Pink	Polystyrene	1000
PBT051	2	Di Potassium EDTA	● Pink	Polystyrene	1000
PBT052	0.5	Lithium Heparin	● Orange	Polystyrene	1000
PBT063	1	Lithium Heparin	● Orange	Polystyrene	1000
PBT0632	2	Lithium Heparin	● Orange	Polystyrene	1000
PBT055	1	Fluoride/Oxalate	● Yellow	Polystyrene	1000
PBT0552	2	Fluoride/Oxalate	● Yellow	Polystyrene	1000
PBT056	2	Sodium Citrate 3.2% (9NC)	● Mauve	Polystyrene	1000