




Capillary Tubes



Glass Capillary Tubes

Pre-Calibrated Micropipettes


- These pipettes are ring marked at a pre-calibrated volume. Volumes are accurate to $\pm 0.25\%$ cv 0.5% .
- The pipettes are colour coded for ease of recognition (ISO 7550). Heparinised version is also available.
- Made from borosilicate glass which has superior chemical properties to avoid contamination of the test medium.



Code	Volume	Colour Code	Na Heparinised	Pack Quantity
GC0029	1-5 μ l	 White		250
GC0063	1-5 μ l	 White	✓	250
GC0031	10 μ l	 Orange		250
GC0064	10 μ l	 Orange	✓	250
GC0065	5-10 μ l	 Orange		250
GC0066	5-10 μ l	 Orange	✓	250
GC0032	20 μ l	 Black		250
GC0067	20 μ l	 Black	✓	250
GC0033	25 μ l	 White / White		250
GC0068	25 μ l	 White / White	✓	250
GC0034	40 μ l	 Red / Red		250
GC0069	40 μ l	 Red / Red	✓	250
GC0056	20-40 μ l	 Red / Red		250
GC0071	20-40 μ l	 Red / Red	✓	250
GC0035	50 μ l	 Green		250
GC0058	50 μ l	 Green	✓	250
GC0057	20-50 μ l	 Green		250
GC0073	20-50 μ l	 Green	✓	250
GC0036	100 μ l	 Blue		250
GC0074	100 μ l	 Blue	✓	250
GC0053	100-200 μ l	 Red		250
GC0047	100-200 μ l	 Red	✓	250

General Purpose Capillaries

- Made from borosilicate glass for superior performance at high temperatures.
- Superior finish and dimensional control makes these products ideal for a range of applications as well as melting point determinations.




Code	Length mm (L)	I.D. mm (d)	O.D. mm	Pack Quantity
GC0044	75	1.0	1.40	1000
GC0042	75	1.2	1.65	500
GC0049	100	1.1	1.55	1000
GC0043	100	1.2	1.65	500
GC0041	100	0.6	1.00	2000
GC0050	125	1.5	1.90	500

* Contact us for details of our extended range, if you do not see the size you require.



Blood Gas Analysis

- Available in 3 sizes.
Heparinised precisely and controlled to a high degree both for size and anticoagulant level.
- Heavy walled with fired tops.
- Stirrers (mixing fleas) and caps (robust and secure for safe handling) available.



Code	Volume	Length mm (L)	Pack Quantity
GC0024	130µl	75	250
GC0025	175µl	100	250
GC0026	220µl	125	250
GC0027	End Cap		500
GC0028	Stirrers (Mixing Fleas)		250

Microhaematocrit Tubes

- Available plain or heparinised (>31.U.) for direct sampling. Taper is less than 0.02mm.
- Contact us for a range of accessories including sealing wax plates.
- Supplied in colour coded dispenser tubes for safe hygienic withdrawal.



Code	Volume	Length mm (L)	Colour Code	Na Heparinised	Pack Quantity
GC0039	75µl	75	Blue		1000 (10 x 100)
GC0040	75µl	75	Red	✓	1000 (10 x 100)
GC0059	40µl	40	Blue		100
GC0060	40µl	40	Red	✓	100

* Please contact us if you do not see the tube you require. Alternative sizes and anti-coagulant types are available to special order.

Accessories available for the above



Code	Description	Pack Quantity
GC0037	Pipette Holder	1
GC0045	Sigillum Wax Plate for Closing Capillary Tubes	10