

Swabs

Culture Swabs

- Wide range of culture swabs available with a choice of materials for the swab shaft.
- CE marked as Class IIa under the Medical Device Directive 93/42/EEC.
- Gamma irradiated.
- Round bottomed tubes for ease of racking.
- Please call if you require a product not shown on these pages, as this only represents a part of the available range.



Code	Description	Case Quantity
SV0015	Culture Swab Plain. Wooden shaft. Single in a tube.	500
SV0018	Culture Swab Plain with Viscose Tip & Polystyrene Shaft. Single in a tube.	500
SV0029	Albu coated Culture Swab Plastic Shaft. Single in a tube.	500
SV0018TS17A	Culture Swab. Viscose Tip. Polypropylene Shaft. Single in a tube.	500
SV0049	Culture Swab with Woodstick. Alginate Tip. Single in a tube.	100
SV0054	Culture Swab. Polystyrene breakable shaft. Polyester tip. Single in a tube	500

Transport Swabs

- Choice of Amies Transport Media and other formulations.
- Colour coded capped swabs with a round bottom tube containing transport media. In a peel pouch.
- CE marked as Class IIa under the Medical Device Directive 93/42/EEC.
- Gamma irradiated.
- Please call if you require a product not shown on these pages, as this only represents a part of the available range.



Code	Description	Qty
SV0021	Swab with Amies Transport Medium. Polystyrene shaft. Viscose tip.	100
SV0023	Swab with Amies Charcoal Transport Medium. Polystyrene shaft. Viscose tip.	100
SV0017T	Mini Tip Swab with Amies Transport Medium. Aluminium shaft. Cotton tip.	100
SV0017TC	Mini Tip Swab with Amies Charcoal Transport Medium. Aluminium shaft. Cotton tip.	100
SV0040	Mini Tip Swab with Amies Charcoal Transport Medium. Plasticised Paper shaft. Cotton Tip.	100

Enviroscreen Swab / Sampling Kit

- A range of environmental and hygiene sampling products.



Code	Description	Quantity
SV0020	5ml Neutralising Buffer. Viscose Tip Swab. Blue polystyrene shaft.	100
SV0048	50mm ² Blue Sponge with Neutralising Buffer in Minigrip Bag	150
SV0041	Dry swab. Blue Polystyrene shaft. Break point at 45mm. Viscose tip. In a tube.	500
SV0033	10ml Neutralising Buffer. Polyester Tip. Polystyrene shaft.	100
SV0063	Pre-Moistened Swab in Neutralising Buffer. Viscose Tip. Blue polystyrene Breakable shaft.	100
SV0064	Liquid MRD Recovery Medium Swab. Viscose Tip.	100
SV0065	Dry Swab. 12" Shaft. Plastic Blue Jumbo Foam Tip. In a tube.	100
SV0066	Dry Swab. Blue Polystyrene Breakable Shaft. Viscose Tip.	1000
SV0072	50cm ² Sponge with Neutralising Buffer with 100ml Plastic Container and Sterile glove.	74
SV0046	Envirostik (2 x 25cm ² Sponge on Stick) with Neutralising Buffer and Medical Grade Bag.	100
SV0068	50cm ² Sponge dosed with Saline (0.9% NaCl) and enclosed Minigrip Sample Bag.	150
SV0070	50cm ² Sponge dosed with Buffered Peptone Water and enclosed Minigrip Sample Bag.	150
SV0071	Swab Kit for the Detection of Pathogenic Listeria. Liquid Media in polypropylene Tube	
	Breakpoint Foam Swab Kit	81

Plating Swabs

- Sterilised, wooden stick plating swabs in a variety of inner pack sizes. For laboratory applications.



Code	Description	Case Quantity
SV0026	Plating Swab. Single Wrapped	1000
SV0016	Plating Swab. Inner pack of 5	1000
SV0024	Plating Swab. Inner pack of 25	1000
SV0014	Plating Swab. Inner pack of 100 (Minigrip Bag)	5000
SV0051	Wooden Applicators 'Orange Sticks'.	1000

Loops and Spreaders

Disposable Microbiological Loops & Spreaders

- Available in bags of 20 in the popular sizes of 1µl, 5µl and 10µl.
- A high precision moulding process produces consistent volumes (certificates on request).
- Mould cavities are highly polished to produce super smooth mouldings.
- Also available are needles (for colony picking) and spreaders.
- Gamma irradiated.
- Other sizes are available on request.



Code	Description	Case Quantity
LOOP01	1µl Loop. Sterile. Inner Sleeve of 20. Clear. Rigid.	1000
LOOP07	1µl Loop. Sterile. Inner Sleeve of 20. Yellow. Soft.	1000
LOOP02	5µl Loop. Sterile. Inner Sleeve of 20. Red. Soft.	1000
LOOP03	10µl Loop. Sterile. Inner Sleeve of 20. Green. Rigid.	1000
LOOP06	10µl Loop. Sterile. Inner Sleeve of 20. White. Soft.	1000
NDL001	Inoculating Needle. Sterile. Inner Sleeve of 20. Rigid.	1000
SPRD01	Spreader. Sterile. Inner Sleeve of 25	1000
SPRD04	Spreader. Sterile. Inner Sleeve of 5	1000

* Individually wrapped and sterilised or alternative packs sizes are available to order.