

Weigh Boats





Anti-static Weigh Boats

Standard

- A full range of disposable weigh boats in white or black.
- Manufactured from semi-rigid plastic, non wettable with smooth surface and rounded corners.
- All boats are flat bottomed with wide rims.
- Inner packs of 250
- Anti-static treatment minimises powder hang-up.
- If you do not see the product that you are looking for then please contact us and we will endeavour to make it for you.



Code	Volume ml	Style	Colour	Case Quantity
WB0004	5ml	Diamond 	Black	1000
WB0003	5ml	Diamond 	White	1000
WB0002	7ml	Square 	Black	1000
WB0001	7ml	Square 	White	1000
WB0032	30ml	Diamond 	Black	500
WB0031	30ml	Diamond 	White	500
WB0030	100ml	Square 	Black	1000
WB0029	100ml	Square 	White	1000
WB0058	100ml	Diamond 	Black	250
WB0057	100ml	Diamond 	White	250
WB0056	250ml	Square 	Black	1000
WB0055	250ml	Square 	White	1000