



Intended use

Specimen receptacle for the containment of samples for laboratory analysis.

Intended users

Professional use only.

Additive

Empty. Option: spoon (for faecal samples).

Handling

Single in vitro use only.

Use for the intended purpose only.

Use appropriate personal protective equipment when handling.

Do not use after the expiry date.

Contamination risk if contents leak – ensure lids and caps are applied correctly.

Do not use if packaging is damaged or product is not sealed.

Only re-cap with correct caps.













Dispose of using appropriate personal protective equipment and act according to applicable regulations.

Tubes maximum centrifugation: 3000 x g for 30 minutes.

Storage

Store in a dry place at room temperatures. Tubes can be refrigerated but it is not necessary to do so.

Keep dry and away from direct sunlight or sources of heat.

 Manufacturer	 Indicates the date after which the product should not be used	 Single use only
 Authorized representative in the European Union	 Manufacturer's batch code	 <i>in vitro</i> diagnostic medical device
 Consult this instruction for use	 Manufacturer's catalogue number	 CE marked (device is in conformity with EU requirements)
 Do not use if package is damaged	 Keep away from sunlight	 Keep dry